



TRANSFER WHEELCHAIRS (DETACHABLE ARMS)

Transit Wheelchair (Detachable Armrests)

Designed for easy patient transfer. Armrests flip up or detach to slide patient from bed to chair.

120 kg

Max User Weight

13 kg

Chair Weight



Product Overview

The ideal choice for caregivers. The movable (detachable/flip-up) armrests allow for 'Side Transfer', enabling the patient to slide smoothly from a bed or sofa onto the chair without lifting. Features 10-inch rear wheels for a compact profile, intended for attendant propulsion.

Ideal Applications

Bed-to-Chair Transfer

Elderly Care

Indoor Transit



Technical Specifications

Feature	Detachable/Flip-up Armrests
Frame	Steel Folding Frame
Max Load (Standard)	120 kg
Rear Wheels	10-inch Solid Rubber-Plastic (Transit)
Seat Width (Standard)	46 cm (Expands to 53 cm)
External Width	59 cm (Fits narrow doors)
Net Weight (Standard)	13 kg
Upholstery	Waterproof Medical Leatherette
Finish	Electrostatic Powder Coating (Rust-Proof)
Braking System	Rear Wheel Parking Brakes
Front Castors	8-inch Rubber

Warranty & After-Sales Support

This product is backed by a manufacturer's warranty of 6 to 12 months. Genuine spare parts, maintenance service, and technical support are provided directly from our Cairo factory. Full tender documentation — datasheets and certificates of origin — is available on request.

Every unit is quality-inspected at our Cairo facility before dispatch.

Contact & Ordering

Phone	+20 109 835 5351	Email	info@aelmasry.com
WhatsApp	+20 109 835 5351	Website	www.aelmasry.com
Alt. Phone	+20 111 114 3969	Address	40 Roshdy Street – Abdeen, Downtown, Cairo, Egypt